



MTA2 Mouse mAb

Catalog No	YP-mAb-19168
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MTA2 MTA1L1 PID
Immunogen	Synthesized peptide derived from human MTA2
Specificity	This antibody detects endogenous levels of MTA2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PID;MTA1L1;A2;MTA2 MTA1L1 PID;Metastasis-associated protein MTA2 (Metastasis-associated 1-like 1) (MTA1-L1 protein) (p53 target protein in deacetylase complex)
Observed Band	75kDa
Calculated Molecular Weight	75kDa
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed.
Function	May be involved in the regulation of gene expression as repressor and activator. The repression might be related to covalent modification of histone proteins.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images