



PSTPIP2 Mouse mAb

Catalog No	YP-mAb-19175
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	PSTPIP2
Immunogen	Synthesized peptide derived from human PSTPIP2
Specificity	This antibody detects endogenous levels of PSTPIP2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MAYP;PSTPIP2;Proline-serine-threonine phosphatase-interacting protein 2 (PEST phosphatase-interacting protein 2)
Observed Band	36 kDa
Calculated Molecular Weight	312 aa, 36 kDa
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	Binds to F-actin. May be involved in regulation of the actin cytoskeleton (By similarity).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images