



INTS10 Mouse mAb

Catalog No	YP-mAb-19184
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	INTS10 C8orf35
Immunogen	Synthesized peptide derived from human INTS10
Specificity	This antibody detects endogenous levels of INTS10 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	INT10;C8orf35;INTS10;INTS10 C8orf35;Integrator complex subunit 10 (Int10)
Observed Band	82kDa
Calculated Molecular Weight	82kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Component of the Integrator (INT) complex, a complex involved in the small nuclear RNAs (snRNA) U1 and U2 transcription and in their 3'-box-dependent processing. The Integrator complex is associated with the C-terminal domain (CTD) of RNA polymerase II largest subunit (POLR2A) and is recruited to the U1 and U2 snRNAs genes (Probable). May be not involved in the recruitment of cytoplasmic dynein to the nuclear envelope by different components of the INT complex .
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images