



WDR92 Mouse mAb

Catalog No	YP-mAb-19210
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	WDR92
Immunogen	Synthesized peptide derived from human WDR92
Specificity	This antibody detects endogenous levels of WDR92 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WD repeat-containing protein Monad;WD repeat-containing protein 92;WD repeat domain 92;Dynein axonemal assembly factor 10;DNAAF10;WDR92;WD repeat-containing protein 92 (WD repeat-containing protein Monad)
Observed Band	40 kDa
Calculated Molecular Weight	357 aa, 40 kDa
Cell Pathway	Dynein axonemal particle .
Tissue Specificity	Widely expressed with the highest expression in testis.
Function	Key assembly factor specifically required for the stability of axonemal dynein heavy chains in cytoplasm.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images