



ICOSL Mouse mAb

Catalog No	YP-mAb-19230
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MON1A SAND1
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	B7h;B7H2;GL50;B7-H2;B7RP1;CD275;ICOSL;LICOS;B7RP-1;ICOS-L
Observed Band	60kDa
Calculated Molecular Weight	33kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables identical protein binding activity. Predicted to be involved in T cell receptor signaling pathway and positive regulation of interleukin-4 production. Located in cytoplasmic ribonucleoprotein granule and plasma membrane
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images