



HILPDA Mouse mAb

Catalog No	YP-mAb-19233
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	
Immunogen	Recombinant protein (or fragment).This information is considered to be commercially sensitive.
Specificity	
Formulation	
Source	Monoclonal,Mouse,IgG
Purification	Affinity purification
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HIG2;HIG-2;C7orf68;HILPDA
Observed Band	Refer to figures
Calculated Molecular Weight	7kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	Enables signaling receptor binding activity. Involved in several processes, including autocrine signaling; cellular response to hypoxia; and positive regulation of lipid storage. Located in several cellular components, including cell surface; lipid droplet; and secretory granule.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images