



CXXC1 Mouse mAb

Catalog No	YP-mAb-19246
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CXXC1 CFP1 CGBP PCCX1 PHF18
Immunogen	Synthesized peptide derived from human CXXC1
Specificity	This antibody detects endogenous levels of CXXC1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CFP1;CGBP;SPP1;PCCX1;PHF18;hCGBP;ZCGPC1;HsT2645;2410002I16Rik;5830420C16Rik;CXXC1;CXXC1 CFP1 CGBP PCCX1 PHF18;CpG-binding protein (CXXC-type zinc finger protein 1) (PHD finger and CXXC domain-containing protein 1)
Observed Band	80 kDa
Calculated Molecular Weight	76 kDa
Cell Pathway	Nucleus speckle . Nucleus . Associated with euchromatin. During mitosis, excluded from condensed chromosomes.
Tissue Specificity	Ubiquitous.
Function	Transcriptional activator that exhibits a unique DNA binding specificity for CpG unmethylated motifs with a preference for CpGG.
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images