



UBE2E1 Mouse mAb

Catalog No	YP-mAb-19288
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	UBE2E1 UBCH6
Immunogen	Synthesized peptide derived from human UBE2E1
Specificity	This antibody detects endogenous levels of UBE2E1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UBCH6;UBE2E1;UBE2E1 UBCH6;Ubiquitin-conjugating enzyme E2 E1 (UbcH6) (Ubiquitin carrier protein E1) (Ubiquitin-protein ligase E1)
Observed Band	21kDa
Calculated Molecular Weight	21kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. Catalyzes the covalent attachment of ISG15 to other proteins. Mediates the selective degradation of short-lived and abnormal proteins. In vitro also catalyzes 'Lys-48'-linked polyubiquitination.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images