



ALG5 Mouse mAb

Catalog No	YP-mAb-19290
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	
Gene Name	ALG5 HSPC149
Immunogen	Synthesized peptide derived from human ALG5
Specificity	This antibody detects endogenous levels of ALG5 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PKD7;bA421P11.2;ALG5;ALG5 HSPC149;Dolichyl-phosphate beta-glucosyltransferase (DoIP-glucosyltransferase) (Asparagine-linked glycosylation protein 5 homolog)
Observed Band	37kDa
Calculated Molecular Weight	37kDa
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type II membrane protein .
Tissue Specificity	Expressed in pancreas, placenta, liver, heart, brain, kidney, skeletal muscle, and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images