



HNMT Mouse mAb

Catalog No	YP-mAb-19322
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HNMT
Immunogen	Synthesized peptide derived from human HNMT
Specificity	This antibody detects endogenous levels of HNMT at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SCN;V19;CP22;CP-22;RI;SRI;Sorcin (22 kDa protein) (CP-22) (V19); HNMT;Histamine N-methyltransferase (HMT) (EC 2.1.1.8)
Observed Band	20kDa
Calculated Molecular Weight	22kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	
Function	Inactivates histamine by N-methylation. Plays an important role in degrading histamine and in regulating the airway response to histamine.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images