



SAMD14 Mouse mAb

Catalog No	YP-mAb-19335
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SAMD14
Immunogen	Synthesized peptide derived from human SAMD14. AA range:297-367
Specificity	This antibody detects endogenous levels of SAMD14 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FZR1 CDH1 FYR FZR KIAA1242;Fizzy-related protein homolog (Fzr) (CDC20-like protein 1) (Cdh1/Hct1 homolog) (hCDH1);SAMD14;Sterile alpha motif domain-containing protein 14 (SAM domain-containing protein 14)
Observed Band	
Calculated Molecular Weight	55kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images