



RPL9 Mouse mAb

Catalog No	YP-mAb-19339
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPL9 OK/SW-cl.103; RPL9P7; RPL9P8; RPL9P9
Immunogen	Synthesized peptide derived from human RPL9. AA range:1-130
Specificity	This antibody detects endogenous levels of RPL9 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SHANK1;SH3 and multiple ankyrin repeat domains protein 1 (Shank1) (Somatostatin receptor-interacting protein) (SSTR-interacting protein) (SSTRIP); RPL9 OK/SW-cl.103;RPL9P7;RPL9P8;RPL9P9;60S ribosomal protein L9
Observed Band	
Calculated Molecular Weight	238kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images