



YPEL1 Mouse mAb

Catalog No	YP-mAb-19340
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	YPEL1 FKSG3
Immunogen	Synthesized peptide derived from human YPEL1. AA range:1-130
Specificity	This antibody detects endogenous levels of YPEL1 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	YPEL1 FKSG3;Protein yippee-like 1
Observed Band	
Calculated Molecular Weight	13kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	May play a role in epithelioid conversion of fibroblasts.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images