



HEBP2 Mouse mAb

Catalog No	YP-mAb-19381
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	HEBP2 C6orf34 SOUL
Immunogen	Synthesized peptide derived from human HEBP2. AA range:72-142
Specificity	This antibody detects endogenous levels of HEBP2 at Human, Mouse, Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HEBP2 C6orf34 SOUL;Heme-binding protein 2 (Placental protein 23) (PP23) (Protein SOUL)
Observed Band	
Calculated Molecular Weight	23kD
Cell Pathway	Cytoplasm . Mitochondrion . Mainly localized to the cytoplasm with a much lower abundance in the mitochondrion. .
Tissue Specificity	Detected in placenta.
Function	Can promote mitochondrial permeability transition and facilitate necrotic cell death under different types of stress conditions.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images