



RNF160 mouse mAb

Catalog No	YP-mAb-19397
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	LTN1 C21orf10 C21orf98 KIAA0714 RNF160 ZNF294 HSPC087
Immunogen	Synthesized peptide derived from human RNF160. AA range:1212-1312
Specificity	This antibody detects endogenous levels of RNF160 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LTN1 C21orf10 C21orf98 KIAA0714 RNF160 ZNF294 HSPC087;E3 ubiquitin-protein ligase listerin (RING finger protein 160) (Zinc finger protein 294)
Observed Band	
Calculated Molecular Weight	194kD
Cell Pathway	Cytoplasm , cytosol .
Tissue Specificity	
Function	E3 ubiquitin-protein ligase. Component of the ribosome quality control complex (RQC) , a ribosome-associated complex that mediates ubiquitination and extraction of incompletely synthesized nascent chains for proteasomal degradation . Ubiquitination leads to VCP/p97 recruitment for extraction and degradation of the incomplete translation product (By similarity) .
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images