



C6orf15 mouse mAb

Catalog No	YP-mAb-19399
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	C6orf15 STG UNQ1840/PRO3566
Immunogen	Synthesized peptide derived from human C6orf15. AA range:6-76
Specificity	This antibody detects endogenous levels of C6orf15 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C6orf15 STG UNQ1840/PRO3566;Uncharacterized protein C6orf15 (Protein STG)
Observed Band	
Calculated Molecular Weight	36kD
Cell Pathway	Secreted , extracellular space , extracellular matrix .
Tissue Specificity	Expressed in skin and tonsils.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images