



DDX10 mouse mAb

Catalog No	YP-mAb-19402
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	DDX10
Immunogen	Synthesized peptide derived from human DDX10. AA range:135-235
Specificity	This antibody detects endogenous levels of DDX10 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DDX10;Probable ATP-dependent RNA helicase DDX10 (DEAD box protein 10)
Observed Band	
Calculated Molecular Weight	96kD
Cell Pathway	
Tissue Specificity	High in testis but widely expressed.
Function	Putative ATP-dependent RNA helicase.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images