



TMEM154 mouse mAb

Catalog No	YP-mAb-19435
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TMEM154
Immunogen	Synthesized peptide derived from human TMEM154. AA range:1-130
Specificity	This antibody detects endogenous levels of TMEM154 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal,Mouse,IgG
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-15°C to -25°C/1 year(Do not lower than -25°C)
Synonyms	TMEM154;Transmembrane protein 154
Observed Band	
Calculated Molecular Weight	20kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

