



CD39 Polyclonal Antibody

Catalog No	YP-Ab-14054
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	ENTPD1
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human ENTDP1. AA range:31-80
Specificity	CD39 Polyclonal Antibody detects endogenous levels of CD39 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ENTPD1;ATP diphosphohydrolase;ATPDase;ATP-DPH;CD39;Ectonucleoside triphosphate diphosphohydrolase 1;NTPDase 1;Ecto-ATP diphosphohydrolase 1;Ecto-ATPDase 1;Ecto-ATPase 1;Ecto-apyrase;Lymphoid cell activation antigen
Observed Band	56kD
Calculated Molecular Weight	56kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed primarily on activated lymphoid cells. Also expressed in endothelial tissues. Isoform 1 and isoform 3 are present in both placenta and umbilical vein, whereas isoform 2 is present in placenta only.
Function	catalytic activity:ATP + 2 H(2)O = AMP + 2 phosphate.,cofactor:Calcium., cofactor:Magnesium.,function:In the nervous system, could hydrolyze ATP and other nucleotides to regulate purinergic neurotransmission. Could also be implicated in the prevention of platelet aggregation. Hydrolyzes ATP and ADP equally well.,PTM:Palmitoylated in the N-terminal part.,PTM:The N-terminus is blocked.,similarity:Belongs to the GDA1/CD39 NTPase family., subunit:Homodimer; disulfide-linked.,tissue specificity:Expressed primarily on activated lymphoid cells. Also expressed in endothelial tissues. The vascular



isoform and the placental isoform II are present in both placenta and umbilical vein, whereas placental isoform I is present in placenta only.,

Background

The protein encoded by this gene is a plasma membrane protein that hydrolyzes extracellular ATP and ADP to AMP. Inhibition of this protein's activity may confer anticancer benefits. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2015],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of HeLa cells using CD39 Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000