



IRAK-4 Monoclonal Antibody

Catalog No	YP-Ab-14163
Isotype	IgG
Reactivity	Human;Mouse;Monkey
Applications	WB;IHC;IF;FCM;ELISA
Gene Name	IRAK4
Immunogen	Purified recombinant fragment of human IRAK-4 expressed in E. Coli.
Specificity	IRAK-4 Monoclonal Antibody detects endogenous levels of IRAK-4 protein.
Formulation	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	WB: 1/500 - 1/2000. IHC: 1/200 - 1/1000. Flow cytometry: 1/200 - 1/400. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1A4C4;EC:2.7.11.1;Interleukin-1 receptor-associated kinase 4;IPD1;IRAK 4; IRAK-4;IRAK4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	extracellular space,intracellular,nucleus,cytoplasm,cytosol,plasma membrane, endosome membrane,
Tissue Specificity	Brain,Spleen,Uterus,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein., cofactor:Magnesium.,disease:Defects in IRAK4 are the cause of IRAK4 deficiency [MIM:607676]. IRAK4 deficiency causes extracellular pyogenic bacterial and fungal infections in otherwise healthy children.,disease:Defects in IRAK4 are the cause of recurrent isolated invasive pneumococcal disease type 1 (IPD1) [MIM:610799]. Recurrent invasive pneumococcal disease (IPD) is defined as two episodes of IPD occurring at least 1 month apart, whether caused by the same or different serotypes or strains. Recurrent IPD occurs in at least 2% of patients in most series, making IPD the most important known risk factor for subsequent IPD.,function:Required for the efficient recruitment of IRAK1 to the IL-1 receptor complex following IL-1 engagement, triggering



intracellular signaling cascades leading to transcriptional up-regulation and mRNA st

Background

This gene encodes a kinase that activates NF-kappaB in both the Toll-like receptor (TLR) and T-cell receptor (TCR) signaling pathways. The protein is essential for most innate immune responses. Mutations in this gene result in IRAK4 deficiency and recurrent invasive pneumococcal disease. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Western Blot analysis using IRAK-4 Monoclonal Antibody against THP-1 (1), HeLa (2), K562 (3), MCF-7 (4), RAW264.7 (5), Jurkat (6) and Cos7 (7) cell lysate.



Immunohistochemistry analysis of paraffin-embedded human lung cancer tissues with DAB staining using IRAK-4 Monoclonal Antibody.



Flow cytometric analysis of Hela cells using IRAK-4 Monoclonal Antibody (blue) and negative control (red).

