



B Raf mouse mAb

Catalog No	YP-Ab-14215
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	braf
Immunogen	Purified recombinant human B Raf protein fragments expressed in E.coli.
Specificity	This antibody detects endogenous levels of B Raf and does not cross-react with related proteins.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	wb 1:1000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Raf-B;braf;FLJ95109;94 kDa B raf protein;B raf 1;B Raf proto oncogene serine threonine protein kinase;B-Raf proto-oncogene serine/threonine-protein kinase (p94);BRAF 1;BRAF_HUMAN;BRAF1;cRmil;MGC126806;MGC138284;Murine sarcoma viral (v-raf) oncogene homolog B1;Murine sarcoma viral v raf oncogene homolog B1;NS7;oncogene BRAF1;p94;Proto-oncogene B-Raf;Proto-oncogene c-Rmil;RAFB 1;RAFB1;RMIL;Serine/threonine-protein kinase B-raf;v raf murine sarcoma viral oncogene homolog B;v raf murine sarcoma viral oncogene homolog B1;v-Raf murine sarcoma viral oncogene homolog B1
Observed Band	85kD
Calculated Molecular Weight	85kD
Cell Pathway	Nucleus . Cytoplasm . Cell membrane . Colocalizes with RGS14 and RAF1 in both the cytoplasm and membranes. .
Tissue Specificity	Brain and testis.
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,cofactor:Binds 2 zinc ions per subunit.,disease:Defects in BRAF are a cause of cardiofaciocutaneous syndrome (CFC syndrome) [MIM:115150]; also known as cardio-facio-cutaneous syndrome. CFC syndrome is characterized by a



distinctive facial appearance, heart defects and mental retardation. Heart defects include pulmonic stenosis, atrial septal defects and hypertrophic cardiomyopathy. Some affected individuals present with ectodermal abnormalities such as sparse, friable hair, hyperkeratotic skin lesions and a generalized ichthyosis-like condition. Typical facial features are similar to Noonan syndrome. They include high forehead with bitemporal constriction, hypoplastic supraorbital ridges, downslanting palpebral fissures, a depressed nasal bridge, and posteriorly angulated ears with prominent helices. The inheritance of CFC synd

Background

This gene encodes a protein belonging to the raf/mil family of serine/threonine protein kinases. This protein plays a role in regulating the MAP kinase/ERKs signaling pathway, which affects cell division, differentiation, and secretion. Mutations in this gene are associated with cardiofaciocutaneous syndrome, a disease characterized by heart defects, mental retardation and a distinctive facial appearance. Mutations in this gene have also been associated with various cancers, including non-Hodgkin lymphoma, colorectal cancer, malignant melanoma, thyroid carcinoma, non-small cell lung carcinoma, and adenocarcinoma of lung. A pseudogene, which is located on chromosome X, has been identified for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot detection of B Raf in MDA-MB-468, T47D, MCF7, HCT116, C2C12, 3T3 and HeLa cell lysates using B Raf mouse mAb (1:1000 diluted). Predicted band size: 87KDa. Observed band size: 87KDa.