



Rsk-1 (phospho Thr359/S363) Polyclonal Antibody

Catalog No	YP-Ab-14451
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	RPS6KA1
Immunogen	The antiserum was produced against synthesized peptide derived from human p90 RSK around the phosphorylation site of Thr359 and Ser363. AA range:331-380
Specificity	Phospho-Rsk-1 (T359/S363) Polyclonal Antibody detects endogenous levels of Rsk-1 protein only when phosphorylated at T359/S363.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/5000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	p RPS6KA1;p RPS6KA1 (Ser380);p RPS6KA1 Ser380;p90RSK (Ser380);Rsk-1;RPS6KA1;Ribosomal protein S6 kinase alpha-1;MAPKAPK1A;RSK1;S6K-alpha-1;90 kDa ribosomal protein S6 kinase 1;p90-RSK 1;p90RSK1;p90S6K;MAP kinase-activated protein kinase 1a;MAPK-activated protein kinase 1a;MAPKAP kinase 1a;MAPKAP
Observed Band	83kD
Calculated Molecular Weight	83kD
Cell Pathway	Nucleus. Cytoplasm.
Tissue Specificity	Colon,Epithelium,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,cofactor:Magnesium.,enzyme regulation:Activated by multiple phosphorylations on threonine and serine residues.,function:Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB.,



PTM:Autophosphorylated on Ser-380, as part of the activation process.,
similarity:Belongs to the protein kinase superfamily.,similarity:Belongs to the
protein kinase superfamily. AGC Ser/Thr protein kinase family. S6 kinase
subfamily.,similarity:Contains 1 AGC-kinase C-terminal domain.,
similarity:Contains 2 protein kinase domains.,subunit:Forms a complex with
either ERK1 or ERK2 in quiescent cells. Transiently dissociates following
mitogenic s

Background

ribosomal protein S6 kinase A1(RPS6KA1) Homo sapiens This gene encodes a
member of the RSK (ribosomal S6 kinase) family of serine/threonine kinases.
This kinase contains 2 nonidentical kinase catalytic domains and phosphorylates
various substrates, including members of the mitogen-activated kinase (MAPK)
signalling pathway. The activity of this protein has been implicated in controlling
cell growth and differentiation. Alternate transcriptional splice variants, encoding
different isoforms, have been characterized. [provided by RefSeq, Jul 2008],

**matters needing
attention**


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For
more information, please consult technical personnel.




Products Images



Enzyme-Linked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phosphopeptide (Phospho-left) and Non-Phosphopeptide (Phospho-right), using p90 RSK (Phospho-Thr359+Ser363) Antibody



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using p90 RSK (Phospho-Thr359+Ser363) Antibody. The picture on the right is blocked with the phospho peptide.



Western blot analysis of lysates from 293 cells treated with PMA 125ng/ml 30', using p90 RSK (Phospho-Thr359+Ser363) Antibody. The lane on the right is blocked with the phospho peptide.