



AK1 Polyclonal Antibody

Catalog No	YP-Ab-14648
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF;ELISA
Gene Name	AK1
Immunogen	The antiserum was produced against synthesized peptide derived from human KAD1 . AA range:101-150
Specificity	AK1 Polyclonal Antibody detects endogenous levels of AK1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	adenylate kinase 1;Adenylate kinase isoenzyme 1;Adenylate monophosphate kinase;AK 1;AK1;ATP-AMP transphosphorylase 1;Myokinase
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm.
Tissue Specificity	Chondrosarcoma Lung Metastasis,Colon,Fetal brain cortex,Retina,Skeletal muscle,Synovial mem
Function	catalytic activity:ATP + AMP = 2 ADP .,disease:Defects in AK1 are the cause of hemolytic anemia due to adenylate kinase deficiency [MIM:612631]., function:Catalyzes the reversible transfer of the terminal phosphate group between ATP and AMP. Small ubiquitous enzyme involved in energy metabolism and nucleotide synthesis that is essential for maintenance and cell growth.,online information:Adenylate kinase entry,polymorphism:This enzyme represents the most common of at least five alleles.,similarity:Belongs to the adenylate kinase family.,subunit:Monomer.,



Background

adenylate kinase 1(AK1) Homo sapiens This gene encodes an adenylate kinase enzyme involved in energy metabolism and homeostasis of cellular adenine nucleotide ratios in different intracellular compartments. This gene is highly expressed in skeletal muscle, brain and erythrocytes. Certain mutations in this gene resulting in a functionally inadequate enzyme are associated with a rare genetic disorder causing nonspherocytic hemolytic anemia. Alternative splicing of this gene results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015],


matters needing attention

Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunofluorescence analysis of HepG2 cells, using KAD1 Antibody . The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using KAD1 Antibody . The picture on the right is blocked with the synthesized peptide.