



HUNK Polyclonal Antibody

Catalog No	YP-Ab-14777
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	HUNK
Immunogen	The antiserum was produced against synthesized peptide derived from human HUNK. AA range:251-300
Specificity	HUNK Polyclonal Antibody detects endogenous levels of HUNK protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	B19;Bstk1;Hormonally up-regulated neu tumor-associated kinase;Hormonally upregulated neu associated kinase;Hormonally upregulated neu tumor associated kinase;HUNK;HUNK_HUMAN;MAKV;Serine/threonine-protein kinase MAK-V
Observed Band	75kD
Calculated Molecular Weight	75kD
Cell Pathway	nucleus,cytoplasm,
Tissue Specificity	Lung adenocarcinoma,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. SNF1 subfamily.,similarity:Contains 1 protein kinase domain.,
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. SNF1 subfamily.,similarity:Contains 1 protein kinase domain.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of the lysates from 293 cells using HUNK antibody.