



HXK I Polyclonal Antibody


Catalog No	YP-Ab-14778
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	HK1
Immunogen	The antiserum was produced against synthesized peptide derived from human HXK1. AA range:31-80
Specificity	HXK I Polyclonal Antibody detects endogenous levels of HXK I protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BB4413;Brain form hexokinase;dea;DrHXK1;EC 2.7.1.1;Glycolytic enzyme; HEXOKIN;hexokinase I;Hexokinase PI;Hexokinase type I;Hexokinase;tumor isozyme;Hexokinase-1;Hexokinase-A;HK I;HK1;HK1 tb;Hk1-s;HK1-ta;HK1-tb; HK1-tc;HKD;HKI;HMSNR;HXK1;HXK1 HUMAN;im:7148527;mHk1-s;wu:fc9d8; wu:fc16e2;wu:fc21e2;wu:fq14b11;zgc:5579;zgc:77618;HXK I
Observed Band	109kD
Calculated Molecular Weight	109kD
Cell Pathway	Mitochondrion outer membrane ; Peripheral membrane protein . Cytoplasm, cytosol . The mitochondrial-binding peptide (MBP) region promotes association with the mitochondrial outer membrane (Probable). Dissociates from the mitochondrial outer membrane following inhibition by N-acetyl-D-glucosamine, leading to relocation to the cytosol (PubMed:27374331). .
Tissue Specificity	Isoform 2: Erythrocyte specific (Ref.6). Isoform 3: Testis-specific (PubMed:10978502). Isoform 4: Testis-specific (PubMed:10978502).




Function	catalytic activity:ATP + D-hexose = ADP + D-hexose 6-phosphate., disease:Defects in HK1 are the cause of hexokinase deficiency [MIM:235700]. Hexokinase deficiency is a rare autosomal recessive disease with nonspherocytic hemolytic anemia as the predominant clinical feature., domain:The N- and C-terminal halves of this hexokinase show extensive sequence similarity to each other. The catalytic activity is associated with the C- terminus while regulatory function is associated with the N-terminus.,enzyme regulation:Hexokinase is an allosteric enzyme inhibited by its product Glc-6-P., miscellaneous:In vertebrates there are four major glucose-phosphorylating isoenzymes, designated hexokinase I, II, III and IV (glucokinase).,online information:Hexokinase entry,pathway:Carbohydrate metabolism; hexose metabolism.,similarity:Belongs to the hexokinase family.,subcellular location:Its hydrophobic N-ter
Background	Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes a ubiquitous form of hexokinase which localizes to the outer membrane of mitochondria. Mutations in this gene have been associated with hemolytic anemia due to hexokinase deficiency. Alternative splicing of this gene results in several transcript variants which encode different isoforms, some of which are tissue- specific. [provided by RefSeq, Apr 2016],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images




Gao, Yunshu, et al. "Overexpression of metabolic markers HK1 and PKM2 contributes to lymphatic metastasis and adverse prognosis in Chinese gastric cancer." *International journal of clinical and experimental pathology* 8.8 (2015): 9264.




Western blot analysis of HXK I Polyclonal Antibody, using Hela, MCF7, SW480, 293T cell, 4° over night, secondary antibody (cat: RS0002) was diluted at 1:10000, 37° 1 hour.




MOUSE-BRAIN



Western Blot analysis of various cells using HXK I
Polyclonal Antibody diluted at 1:1000



Western Blot analysis of HeLa cells using HXK I
Polyclonal Antibody diluted at 1:1000



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody was pre-absorbed by immunogen peptide.