



PKD1/2/3 Polyclonal Antibody

Catalog No	YP-Ab-14940
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA;IHC
Gene Name	KPCD1
Immunogen	The antiserum was produced against synthesized peptide derived from human PKD1/2/3/PKC mu. AA range:706-755
Specificity	PKD1/2/3 Polyclonal Antibody detects endogenous levels of PKD1/2/3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZp586E82;HSPC187;KPCD2_HUMAN;nPKC D2;nPKC-D2;PKD2;PRKD 2;Prkd2;Protein kinase D2;Serine/threonine protein kinase D2;Serine/threonine-protein kinase D2;KPCD1_HUMAN;nPKC D1;nPKC mu;nPKC-D1;nPKC-mu;nPKCD1;nPKCmu;PKC;PKC MU;PKCM;PKCmu;PKD 1;PKD;PKD1;PRKCM;PRKD 1;Prkd1;Protein kinase C mu;Protein kinase C mu type;Protein kinase D;Protein kinase D1;Serine/threonine protein kinase D1;Serine/threonine-protein kinase D1;EPK 2;EPK2;KPCD3_HUMAN;nPKC nu;nPKC-nu;nPKCnu;nu;PKCnu;PKD 3;PKD3;PRK D3;PRKCN;PRKD 3;Prkd3;Protein kinase C;Protein kinase C nu;Protein kinase C nu type;Protein kinase D3;Protein kinase EPK 2;Protein kinase EPK2;Serine threonine protein kinase;Serine/threonine protein kinase D3;Serine/threonine-protein kinase D3;PKD1/2/3;KPCD1;Serine/threonine-protein kinas
Observed Band	115kD
Calculated Molecular Weight	115kD
Cell Pathway	Cytoplasm . Cell membrane . Golgi apparatus, trans-Golgi network . Translocation to the cell membrane is required for kinase activation.
Tissue Specificity	Placenta,Testis,



Function

catalytic activity:ATP + a protein = ADP + a phosphoprotein.,enzyme regulation:Activated by diacylglycerol and phorbol esters.,function:Calcium-independent, phospholipid-dependent, serine- and threonine-specific kinase involved in resistance to oxidative stress.,PTM:Phosphorylation of Ser-738 and/or Ser-742 in activated PKD is mediated by transphosphorylation (By similarity). Phosphorylation of Tyr-463 mediated by the Src/Abl pathway in response to oxidative stress activates the kinase.,similarity:Belongs to the protein kinase superfamily.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. PKD subfamily.,similarity:Contains 1 PH domain.,similarity:Contains 1 protein kinase domain.,similarity:Contains 2 phorbol-ester/DAG-type zinc fingers.,subunit:Interacts (via N-terminus) with ADAP1/CENTA1. Interacts with Src.,

Background

PRKD1 is a serine/threonine kinase that regulates a variety of cellular functions, including membrane receptor signaling, transport at the Golgi, protection from oxidative stress at the mitochondria, gene transcription, and regulation of cell shape, motility, and adhesion (summary by Eiseler et al., 2009 [PubMed 19329994]).[supplied by OMIM, Nov 2010],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from Jurkat cells, using PKD1/2/3/PKC mu Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).