



PPP4R1L Polyclonal Antibody


Catalog No	YP-Ab-14957
Isotype	IgG
Reactivity	Human;Monkey
Applications	WB;ELISA
Gene Name	PPP4R1L
Immunogen	The antiserum was produced against synthesized peptide derived from human PPP4R1L. AA range:191-240
Specificity	PPP4R1L Polyclonal Antibody detects endogenous levels of PPP4R1L protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C2orf192;PRO185;Protein phosphatase 4 regulatory subunit 1 like;PPP4R1L; Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like; C20orf192
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	
Tissue Specificity	Fetal liver,Thymus,
Function	function:May be a regulatory subunit of serine/threonine-protein phosphatase 4., similarity:Contains 4 HEAT repeats.,
Background	function:May be a regulatory subunit of serine/threonine-protein phosphatase 4., similarity:Contains 4 HEAT repeats.,
matters needing attention	Avoid repeated freezing and thawing!




Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.


Products Images



Western Blot analysis of various cells using PPP4R1L Polyclonal Antibody



Western blot analysis of lysates from COLO, HUVEC, HepG2, Jurkat, and COS cells, using PPP4R1L Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using PPP4R1L antibody.