



TNNI3K Polyclonal Antibody

Catalog No	YP-Ab-15021
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	TNNI3K
Immunogen	The antiserum was produced against synthesized peptide derived from human TNNI3K. AA range:301-350
Specificity	TNNI3K Polyclonal Antibody detects endogenous levels of TNNI3K protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CARK;CCDD;TNNI3K;Serine/threonine-protein kinase TNNI3K;Cardiac ankyrin repeat kinase;Cardiac troponin I-interacting kinase;TNNI3-interacting kinase
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cytoplasm . Expressed at lower levels in the cytoplasm.
Tissue Specificity	Highly expressed in both adult and fetal heart.
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,catalytic activity:GTP + beta-L-fucose 1-phosphate = diphosphate + GDP-L-fucose., cofactor:Magnesium.,function:Catalyzes the formation of GDP-L-fucose from GTP and L-fucose-1-phosphate. Functions as a salvage pathway to reutilize L-fucose arising from the turnover of glycoproteins and glycolipids.,function:May play a role in cardiac physiology.,PTM:Autophosphorylated.,similarity:Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family. MAP kinase kinase kinase subfamily.,similarity:Contains 1 protein kinase domain., similarity:Contains 10 ANK repeats.,subcellular location:Expressed at lower levels in the cytoplasm.,subunit:Interacts with TNNI3, ACTC, ACTA1, MYBPC3, AIP, BABP3 and HADHB.,tissue specificity:Expressed in many tissues.,tissue



specificity: Highly expressed in both adult and fetal heart.,

Background

This gene encodes a protein that belongs to the MAP kinase kinase kinase (MAPKKK) family of protein kinases. The protein contains ankyrin repeat, protein kinase and serine-rich domains and is thought to play a role in cardiac physiology. [provided by RefSeq, Sep 2012],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Mouse-kidney



Western blot analysis of Mouse-kidney lysis using TNNI3K antibody. Antibody was diluted at 1:500. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using TNNI3K Antibody. The picture on the right is blocked with the synthesized peptide.