



AMPK β 2 rabbit pAb

Catalog No	YP-Ab-15056
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB; ELISA
Gene Name	PRKAB2
Immunogen	Synthesized peptide derived from human AMPK β 2 AA range: 130-210
Specificity	This antibody detects endogenous levels of Human AMPK β 2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	PRKAB2;5'-AMP-activated protein kinase subunit beta-2;AMPK;AMPK subunit beta 2;AMPK subunit beta-2;AMPK β 2;5'-AMP-activated protein kinase subunit beta-2 (AMPK subunit beta-2)
Observed Band	34 kDa
Calculated Molecular Weight	30 kDa
Cell Pathway	
Tissue Specificity	
Function	fatty acid metabolic process, fatty acid biosynthetic process, lipid biosynthetic process, regulation of cellular ketone metabolic process, organic acid biosynthetic process, regulation of lipid metabolic process, regulation of fatty acid metabolic process, regulation of fatty acid oxidation, carboxylic acid biosynthetic process,
Background	function:AMPK is responsible for the regulation of fatty acid synthesis by phosphorylation of acetyl-CoA carboxylase. Also regulates cholesterol synthesis via phosphorylation and inactivation of hydroxymethylglutaryl-CoA reductase and hormone-sensitive lipase. This is a regulatory subunit, may be a positive regulator of AMPK activity. It may also serve as an adapter molecule for the



catalytic alpha-subunit.,PTM:Phosphorylated when associated with the catalytic subunit.,similarity:Belongs to the 5'-AMP-activated protein kinase beta subunit family.,subunit:Heterotrimer of an alpha catalytic subunit, a beta and a gamma non-catalytic regulatory subunits.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images