



MutS Protein Homolog 6(MSH6) (PT0171) mouse mAb

Catalog No	YP-Ab-15278
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF/ICC
Gene Name	MSH6 GTBP
Immunogen	Synthesized peptide derived from human MutS Protein Homolog 6(MSH6)
Specificity	This antibody detects endogenous levels of human MutS Protein Homolog 6(MSH6). Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse/IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, WB 1:500-2000, IF 1:500-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MSH6 GTBP;DNA mismatch repair protein Msh6 (hMSH6) (G/T mismatch-binding protein) (GTBP) (GTMBP) (MutS-alpha 160 kDa subunit) (p160)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome . Associates with H3K36me3 via its PWWP domain.
Tissue Specificity	Epithelium,Placenta,Pooled,Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images