



Mesothelin mouse mAb(ABT-MSLN)

Catalog No	YP-Ab-15470
Isotype	IgG
Reactivity	Human
Applications	IHC;IF
Gene Name	MSLN MPF
Immunogen	Synthesized peptide derived from human Mesothelin
Specificity	This antibody detects endogenous levels of human Mesothelin . Heat-induced epitope retrieval (HIER) Citrate buffer of pH6.0 was highly recommended as antigen repair method in paraffin section
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Mouse, Monoclonal/IgG2b, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MSLN;MPF Mesothelin;cleaved form;MPF;Mesothelin;Megakaryocyte-potentiating factor;MSLN MPF;Mesothelin (CAK1 antigen;Pre-pro-megakaryocyte-potentiating factor) [Cleaved into: Megakaryocyte-potentiating factor (MPF);cleaved form]
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane; Lipid-anchor, GPI-anchor. Golgi apparatus.; [Megakaryocyte-potentiating factor]: Secreted.; [Isoform 3]: Secreted.
Tissue Specificity	Expressed in lung. Expressed at low levels in heart, placenta and kidney. Expressed in mesothelial cells. Highly expressed in mesotheliomas, ovarian cancers, and some squamous cell carcinomas (at protein level).
Function	disease:Antibodies against MSLN are detected in patients with mesothelioma and ovarian cancer.;function:Megakaryocyte-potentiating factor (MPF) potentiates megakaryocyte colony formation in vitro.;function:Membrane-anchored forms may play a role in cellular adhesion.;PTM:Both MPF and the cleaved form of mesothelin are N-glycosylated.;PTM:Proteolytically cleaved by a furin-like convertase to generate megakaryocyte-potentiating factor (MPF), and



the cleaved form of mesothelin.,similarity:Belongs to the mesothelin family.,
subunit:Interacts with MUC16.,tissue specificity:Expressed in lung. Expressed at
low levels in heart, placenta and kidney. Expressed in mesothelial cells. Highly
expressed in mesotheliomas, ovarian cancers, and some squamous cell
carcinomas (at protein level).,

Background

This gene encodes a preproprotein that is proteolytically processed to generate two protein products, megakaryocyte potentiating factor and mesothelin. Megakaryocyte potentiating factor functions as a cytokine that can stimulate colony formation of bone marrow megakaryocytes. Mesothelin is a glycosylphosphatidylinositol-anchored cell-surface protein that may function as a cell adhesion protein. This protein is overexpressed in epithelial mesotheliomas, ovarian cancers and in specific squamous cell carcinomas. Alternative splicing results in multiple transcript variants, at least one of which encodes an isoform that is proteolytically processed. [provided by RefSeq, Feb 2016],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images



Immunohistochemical analysis of paraffin-embedded Mesothelioma. 1, Antibody was diluted at 1:200(4° overnight). 2, Citrate buffer of pH6.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 30min).



Immunohistochemical analysis of paraffin-embedded Ovarian Serous Carcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Citrate buffer of pH6.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 30min).



Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Mesothelin (ABT-MSLN) antibody. The HRP-conjugated Goat anti-Mouse IgG(H + L) antibody was used to detect the antibody. Lane 1: HeLa Predicted band size: 69kDa Observed band size: 69kDa