



MUC5AC mouse mAb(ABT350)

Catalog No	YP-Ab-15604
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	MUC5AC MUC5
Immunogen	Synthesized peptide derived from human MUC5AC
Specificity	This antibody detects endogenous levels of MUC5AC at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.161% sodium azide.
Source	Monoclonal, Mouse/IgG2a, kappa
Purification	he antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	0.56mg/mL
Purity	≥90%
Storage Stability	-20°C/1 year
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in larger amounts in spleen and lung. Detected in liver, heart, small intestine, colon, thymus, prostate, skeletal muscle, peripheral blood leukocytes, kidney and testis. No expression in the brain.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images