



BSX Polyclonal Antibody

Catalog No	YP-Ab-15750
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	BSX
Immunogen	Synthesized peptide derived from BSX . at AA range: 150-230
Specificity	BSX Polyclonal Antibody detects endogenous levels of BSX protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Brain specific homeobox;Brain specific homeobox protein homolog;mouse; homolog of;BSH;BSX1;BSX;Brain-specific homeobox protein homolog
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:DNA binding protein that function as transcriptional activator. Is essential for normal postnatal growth and nursing. Is an essential factor for neuronal neuropeptide Y and agouti-related peptide function and locomotory behavior in the control of energy balance.,online information:Of fidgets and food - Issue 85 of August 2007,similarity:Belongs to the distal-less homeobox family., similarity:Contains 1 homeobox DNA-binding domain.,
Background	function:DNA binding protein that function as transcriptional activator. Is essential for normal postnatal growth and nursing. Is an essential factor for neuronal neuropeptide Y and agouti-related peptide function and locomotory behavior in the control of energy balance.,online information:Of fidgets and food - Issue 85 of August 2007,similarity:Belongs to the distal-less homeobox family.,

similarity:Contains 1 homeobox DNA-binding domain.,


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of 293T cells using BSX Polyclonal Antibody diluted at 1:500