



Oct-6 Polyclonal Antibody

Catalog No	YP-Ab-15790
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	POU3F1
Immunogen	The antiserum was produced against synthesized peptide derived from human OCT6. AA range:311-360
Specificity	Oct-6 Polyclonal Antibody detects endogenous levels of 40822 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OCT6;OTF6;SCIP;POU3F1;45205;POU domain class 3 transcription factor 1; POU domain;class 3;transcription factor 1;Octamer-binding protein 6;Oct-6; Octamer-binding transcription factor 6;OTF-6;POU domain transcription factor SCIP
Observed Band	45kD
Calculated Molecular Weight	45kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in embryonal stem cells and in the developing brain.
Function	function:Transcription factor that binds to the octamer motif (5'-ATTGTCAT-3'). Thought to be involved in early embryogenesis and neurogenesis., similarity:Belongs to the POU transcription factor family. Class-3 subfamily., similarity:Contains 1 homeobox DNA-binding domain.,similarity:Contains 1 POU-specific domain.,tissue specificity:Expressed in embryonal stem cells and in the developing brain.,



Background

function: Transcription factor that binds to the octamer motif (5'-ATTTGCAT-3'). Thought to be involved in early embryogenesis and neurogenesis., similarity: Belongs to the POU transcription factor family. Class-3 subfamily., similarity: Contains 1 homeobox DNA-binding domain., similarity: Contains 1 POU-specific domain., tissue specificity: Expressed in embryonal stem cells and in the developing brain.,


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of various cells using Oct-6 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).



Western blot analysis of lysates from A549 cells, using OCT6 Antibody. The lane on the right is blocked with the synthesized peptide.