




ITIH4 (80k, Cleaved-Arg689) rabbit pAb

Catalog No	YP-Ab-15818
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC
Gene Name	ITIH4 IHRP ITIHL1 PK120 PRO1851
Immunogen	Synthesized peptide derived from human ITIH4 (80k, Cleaved-Arg689)
Specificity	This antibody detects endogenous levels of Human ITIH4 (80k, Cleaved-Arg689)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	35 kDa inter-alpha-trypsin inhibitor heavy chain H4;70 kDa inter-alpha-trypsin inhibitor heavy chain H4;IHRP;Inter-alpha-inhibitor heavy chain 4;Inter-alpha-trypsin inhibitor family heavy chain-related protein;ITIH4;ITIH4 IHRP ITIHL1 PK120 PRO1851;ITIH4 (80k;Cleaved-Arg689);Inter-alpha-trypsin inhibitor heavy chain H4 (ITI heavy chain H4;ITI-HC4;Plasma kallikrein sensitive glycoprotein 120;Gp120;PK-120) [Cleaved into: 70 kDa inter-alpha-trypsin inhibitor heavy chain H4;35 kDa inter-alpha-trypsin inhibitor heavy chain H4]
Observed Band	80 100kD
Calculated Molecular Weight	80 100kD
Cell Pathway	Secreted.
Tissue Specificity	Liver specific.
Function	ERROR 23
Background	ERROR 23
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).