



EFNA3 (Cleaved-Gly214) rabbit pAb

Catalog No	YP-Ab-15820
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	EFNA3 EFL2 EPLG3 LERK3
Immunogen	Synthesized peptide derived from human EFNA3 (Cleaved-Gly214)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat EFNA3 (Cleaved-Gly214, protein was cleaved amino acid sequence between 214-215)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EFL2;EFNA3;EPLG3;LERK3;EFL 2;Ephrin-A3;EFNA3 EFL2 EPLG3 LERK3; EFNA3 (Cleaved-Gly214);Ephrin-A3 (EFL-2;EHK1 ligand;EHK1-L;EPH-related receptor tyrosine kinase ligand 3;LERK-3)
Observed Band	26-35 kDa
Calculated Molecular Weight	238 aa, 26 kDa
Cell Pathway	Cell membrane; Lipid-anchor, GPI-anchor.
Tissue Specificity	Expressed in brain, skeletal muscle, spleen, thymus, prostate, testis, ovary, small intestine, and peripheral blood leukocytes.
Function	cell-cell signaling,
Background	similarity:Belongs to the ephrin family.,tissue specificity:Expressed in brain, skeletal muscle, spleen, thymus, prostate, testis, ovary, small intestine, and peripheral blood leukocytes.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images