




CTRP8 Polyclonal Antibody

Catalog No	YP-Ab-15894
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	C1QTNF8
Immunogen	Synthesized peptide derived from the C-terminal region of human CTRP8.
Specificity	CTRP8 Polyclonal Antibody detects endogenous levels of CTRP8 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1q and tumor necrosis factor related protein 8;C1q/TNF-related protein 8; C1QT8_HUMAN;C1QTNF8;Complement C1q tumor necrosis factor-related protein 8;CTRP8;UNQ5829;UNQ5829/PRO19648
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed predominantly in lung and testis (PubMed:19666007). Expressed in astrocytes (PubMed:24014093).
Function	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,
Background	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).