



GDF-15 Polyclonal Antibody

Catalog No	YP-Ab-15914
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	GDF15
Immunogen	The antiserum was produced against synthesized peptide derived from human GDF15. AA range:31-80
Specificity	GDF-15 Polyclonal Antibody detects endogenous levels of GDF-15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GDF15;GDF 15;growth differentiation factor 15;Growth/differentiation factor 15; MIC-1;GDF-15;MIC1;PDF;PLAB;PTGFB;Macrophage inhibitory cytokine 1; NSAID-activated gene 1 protein;NAG-1;NSAID-regulated gene 1 protein;NRG-1; Placental TGF-beta;Placental bone morphogenetic
Observed Band	37kD
Calculated Molecular Weight	37kD
Cell Pathway	Secreted .
Tissue Specificity	Highly expressed in placenta, with lower levels in prostate and colon and some expression in kidney (PubMed:9348093). Detected in plasma (at protein level) (PubMed:28572090, PubMed:29046435).
Function	similarity:Belongs to the TGF-beta family.,subunit:Homodimer; disulfide-linked., tissue specificity:Highly expressed in placenta, with lower levels in prostate and colon and some expression in kidney.,
Background	growth differentiation factor 15(GDF15) Homo sapiens This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to



recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. The protein is expressed in a broad range of cell types, acts as a pleiotropic cytokine and is involved in the stress response program of cells after cellular injury. Increased protein levels are associated with disease states such as tissue hypoxia, inflammation, acute injury and oxidative stress. [provided by RefSeq, Aug 2016],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis of various cells using GDF-15 Polyclonal Antibody diluted at 1:1000



Western blot analysis of lysates from 293 and HepG2 cells, using GDF15 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using GDF15 antibody.

