



GDF-9 Polyclonal Antibody

Catalog No	YP-Ab-15916
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;IHC;IF;ELISA
Gene Name	GDF9
Immunogen	The antiserum was produced against synthesized peptide derived from human GDF-9. AA range:273-322
Specificity	GDF-9 Polyclonal Antibody detects endogenous levels of GDF-9 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FL;Flt 3 ligand;Flt 3L;Flt3 L;FLT3 LG;Flt3 ligand;Flt3L;FLT3L_HUMAN;Flt3lg; Fms related tyrosine kinase 3 ligand;Fms-related tyrosine kinase 3 ligand;SL cytokine;GDF-9;GDF9;Growth/differentiation factor 9
Observed Band	57kD
Calculated Molecular Weight	57kD
Cell Pathway	Secreted .
Tissue Specificity	Expressed in ovarian granulosa cells. Present in oocytes of primary follicles (at protein level).
Function	function:Required for ovarian folliculogenesis.,similarity:Belongs to the TGF-beta family.,subunit:Homodimer or heterodimer (Potential). But, in contrast to other members of this family, cannot be disulfide-linked.,
Background	This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked



homodimer. This protein regulates ovarian function. Reduced expression of this gene may be associated with polycystic ovary syndrome and mutations in this gene may be more common in mothers of dizygotic twins. [provided by RefSeq, Jul 2016],

matters needing attention


Avoid repeated freezing and thawing!

Usage suggestions


This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.




Products Images




Western Blot analysis of various cells using GDF-9 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded Human lung cancer. Antibody was diluted at 1:100 (4° overnight). High-pressure and temperature Tris-EDTA, pH 8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.



Immunohistochemistry analysis of GDF-9 antibody in paraffin-embedded human liver carcinoma tissue.



Western blot analysis of lysate from COLO205 cells, using GDF-9 antibody.