



IGF-I Polyclonal Antibody

Catalog No	YP-Ab-15927
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	IHC;IF;ELISA
Gene Name	IGF1
Immunogen	The antiserum was produced against synthesized peptide derived from human IGF-I. AA range:100-149
Specificity	IGF-I Polyclonal Antibody detects endogenous levels of IGF-I protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/10000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IBP1;IGF I;IGF1;IGF1A;IGF1B;IGFI;Insulin like growth factor I;Mechano growth factor;MGF;Somatomedin C;IGF-I;Insulin-like growth factor I;IGF-IA;insulin-like growth factor 1;insulin-like growth factor IA
Observed Band	
Calculated Molecular Weight	
Cell Pathway	extracellular region,extracellular space,plasma membrane,insulin-like growth factor binding protein complex,platelet alpha granule lumen,alphav-beta3 integrin-IGF-1-IGF1R complex,insulin-like growth factor ternary complex, exocytic vesicle,
Tissue Specificity	Brain,Liver,PCR rescued clones,
Function	disease:Defects in IGF1 are the cause of insulin-like growth factor I deficiency (IGF1 deficiency) [MIM:608747]. IGF1 deficiency is an autosomal recessive disorder characterized by growth retardation, sensorineural deafness and mental retardation.,function:The insulin-like growth factors, isolated from plasma, are structurally and functionally related to insulin but have a much higher growth-promoting activity.,online information:Insulin-like growth factor 1 entry,online information:The Singapore human mutation and polymorphism database,



similarity:Belongs to the insulin family.,

Background

The protein encoded by this gene is similar to insulin in function and structure and is a member of a family of proteins involved in mediating growth and development. The encoded protein is processed from a precursor, bound by a specific receptor, and secreted. Defects in this gene are a cause of insulin-like growth factor I deficiency. Alternative splicing results in multiple transcript variants encoding different isoforms that may undergo similar processing to generate mature protein. [provided by RefSeq, Sep 2015],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded Human pancreas. Antibody was diluted at 1:100(4° overnight). High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody was pre-absorbed by immunogen peptide.



Immunohistochemistry analysis of IGF-I antibody in paraffin-embedded human pancreas tissue.