



Osteoglycin Polyclonal Antibody

Catalog No	YP-Ab-16011
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	OGN
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human OGN. AA range:231-280
Specificity	Osteoglycin Polyclonal Antibody detects endogenous levels of Osteoglycin protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OGN;SLRR3A;Osteoinductive factor;OIF;OG;Osteoglycin;Mimecan
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Bone.
Function	function:Induces bone formation in conjunction with TGF-beta-1 or TGF-beta-2., similarity:Belongs to the small leucine-rich proteoglycan (SLRP) family. Class III subfamily.,similarity:Contains 7 LRR (leucine-rich) repeats.,tissue specificity:Bone.,
Background	This gene encodes a member of the small leucine-rich proteoglycan (SLRP) family of proteins. The encoded protein induces ectopic bone formation in conjunction with transforming growth factor beta and may regulate osteoblast differentiation. High expression of the encoded protein may be associated with elevated heart left ventricular mass. Alternative splicing results in multiple



transcript variants. [provided by RefSeq, Jan 2016],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis of L929 cells using Osteoglycin Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



Western blot analysis of lysate from L929 cells, using OGN Antibody.