



HRP-3 Polyclonal Antibody


Catalog No	YP-Ab-16040
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC;IF;ELISA
Gene Name	HDGFRP3
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human HDGFRP3. AA range:131-180
Specificity	HRP-3 Polyclonal Antibody detects endogenous levels of HRP-3 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HDGF2;HRP-3;HDGF-2;CGI-142;HDGFRP3;Hepatoma-derived growth factor-related protein 3;Hepatoma-derived growth factor 2
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	Nucleus .
Tissue Specificity	Detected in testis, heart, spinal cord and brain.
Function	function:Enhances DNA synthesis and may play a role in cell proliferation., similarity:Belongs to the HDGF family.,similarity:Contains 1 PWWP domain., tissue specificity:Detected in testis, heart, spinal cord and brain.,
Background	function:Enhances DNA synthesis and may play a role in cell proliferation., similarity:Belongs to the HDGF family.,similarity:Contains 1 PWWP domain., tissue specificity:Detected in testis, heart, spinal cord and brain.,
matters needing attention	Avoid repeated freezing and thawing!




Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.


Products Images



Western Blot analysis of NIH-3T3, PC12 cells using HRP-3 Polyclonal Antibody. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded human-Lymph-nodes, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-Lymph-nodes, antibody was diluted at 1:100