




RASSF1 Monoclonal Antibody

Catalog No	YP-Ab-16105
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	RASSF1
Immunogen	Purified recombinant human RASSF1 protein fragments expressed in E.coli.
Specificity	RASSF1 Monoclonal Antibody detects endogenous levels of RASSF1 protein.
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/1000 - 1/2000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	123F2;Cardiac specific ras association domain family 1 protein;cardiac-specific ras association domain family 1 protein;NORE2A;OTTHUMP19553; OTTHUMP19554;OTTHUMP19555;OTTHUMP19557;Pancreas specific ras association domain family 1 protein;Ras association (RalGDS/AF 6) domain family member 1;Ras association domain-containing protein 1;RASf1_HUMAN; Rassf1;RASSF1A;RDA32;REH3P21;Tumor suppressor protein RDA32; WUGSC:H_LUCA12.5
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform A]: Cytoplasm, cytoskeleton. Cytoplasm, cytoskeleton, microtubule organizing center, centrosome. Cytoplasm, cytoskeleton, spindle. Cytoplasm, cytoskeleton, spindle pole. Nucleus. Localizes to cytoplasmic microtubules during interphase, to bipolar centrosomes associated with microtubules during prophase, to spindle fibers and spindle poles at metaphase and anaphase, to the midzone during early telophase, and to the midbody in late telophase in cells. Colocalizes with MDM2 in the nucleus.; [Isoform C]: Nucleus. Predominantly nuclear.



Tissue Specificity	Isoform A and isoform C are ubiquitously expressed in all tissues tested, however isoform A is absent in many corresponding cancer cell lines. Isoform B is mainly expressed in hematopoietic cells.
Function	function: Potential tumor suppressor. Required for death receptor-dependent apoptosis. Mediates activation of STK4 during Fas-induced apoptosis. When associated with MOAP1, promotes BAX conformational change and translocation to mitochondrial membranes in response to TNF and TNFSF10 stimulation. Isoform A interacts with CDC20, an activator of the anaphase-promoting complex, APC, resulting in the inhibition of APC activity and mitotic progression. Inhibits proliferation by negatively regulating cell cycle progression at the level of G1/S-phase transition by regulating accumulation of cyclin D1 protein. Isoform C has been shown not to perform these roles, no function has been identified for this isoform.,similarity: Contains 1 phorbol-ester/DAG-type zinc finger.,similarity: Contains 1 Ras-associating domain.,similarity: Contains 1 SARAH domain.,subcellular location: Localizes to cytoplasmic mic
Background	This gene encodes a protein similar to the RAS effector proteins. Loss or altered expression of this gene has been associated with the pathogenesis of a variety of cancers, which suggests the tumor suppressor function of this gene. The inactivation of this gene was found to be correlated with the hypermethylation of its CpG-island promoter region. The encoded protein was found to interact with DNA repair protein XPA. The protein was also shown to inhibit the accumulation of cyclin D1, and thus induce cell cycle arrest. Several alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported. [provided by RefSeq, May 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western Blot analysis using RASSF1 Monoclonal Antibody against mouse liver, mouse brain lysate.