



ARHGEF16 Polyclonal Antibody

Catalog No	YP-Ab-16139
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ARHGEF16
Immunogen	Synthesized peptide derived from ARHGEF16 . at AA range: 140-220
Specificity	ARHGEF16 Polyclonal Antibody detects endogenous levels of ARHGEF16 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NBR;GEF16;Ephexin-4;EPHEXIN4;Ephexin 4;ARHGEF16;Rho guanine nucleotide exchange factor 16
Observed Band	80 kDa
Calculated Molecular Weight	80 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Lung,PNS,
Function	similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,similarity:Contains 1 SH3 domain.,
Background	Although the specific function of this protein is not known yet, it is thought to be involved in protein-protein and protein-lipid interactions. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images