




GBP5 Polyclonal Antibody

Catalog No	YP-Ab-16162
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;WB;ELISA
Gene Name	GBP5
Immunogen	Synthesized peptide derived from the Internal region of human GBP5.
Specificity	GBP5 Polyclonal Antibody detects endogenous levels of GBP5 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 IHC: 1/100 - 1/300. ELISA: 1/20000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.6.5.-;GBP-5;GBP-TA antigen;GTP-binding protein 5;Guanine nucleotide-binding protein 5;GBP5;Guanylate-binding protein 5
Observed Band	66kD
Calculated Molecular Weight	66kD
Cell Pathway	Cytoplasm . Golgi apparatus membrane ; Lipid-anchor . Cytoplasmic vesicle .
Tissue Specificity	Expressed in peripheral blood monocytes (at protein level).
Function	similarity:Belongs to the GBP family.,
Background	similarity:Belongs to the GBP family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human small intestinal carcinoma tissue. 1, primary Antibody was diluted at 1:200 (4° overnight). 2, Sodium citrate pH 6.0 was used for antigen retrieval (>98° C, 20min). 3, Secondary antibody was diluted at 1:200