




Rab 33A Polyclonal Antibody


Catalog No	YP-Ab-16197
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA;IHC
Gene Name	RAB33A
Immunogen	The antiserum was produced against synthesized peptide derived from human RAB33A. AA range:132-181
Specificity	Rab 33A Polyclonal Antibody detects endogenous levels of Rab 33A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Small GTP-binding protein S10;Small GTP binding protein S10;Ras-related protein Rab-33A;Ras related protein Rab 33A;RABS10;Rab 33A;RAB33A
Observed Band	27kD
Calculated Molecular Weight	27kD
Cell Pathway	Cell membrane ; Lipid-anchor ; Cytoplasmic side .
Tissue Specificity	Expressed only in lymphoid cell lines.
Function	similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Expressed only in lymphoid cell lines.,
Background	The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. It is GTP-binding protein and may be involved in vesicle transport. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH 9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45 min).