



Rab 6C Polyclonal Antibody

Catalog No	YP-Ab-16206
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;IHC
Gene Name	RAB6C
Immunogen	The antiserum was produced against synthesized peptide derived from human RAB6C. AA range:165-214
Specificity	Rab 6C Polyclonal Antibody detects endogenous levels of Rab 6C protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RAB 6C;Rab6 like protein WTH3;Rab6-like protein WTH3;RAB6C;RAB6C member RAS oncogene family;RAB6C_HUMAN;Ras related protein Rab 6C; Ras-related protein Rab-6C;WTH3
Observed Band	28kD
Calculated Molecular Weight	28kD
Cell Pathway	Nucleus . Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome .
Tissue Specificity	Highest levels are found in fetal and adult brain, prostate, testis and spinal cord. Undetectable expression in adrenal gland, skeletal muscle, bone marrow, fetal, and adult liver, heart, salivary gland, and trachea. Detected in the HEK293, HEK293T, LNCaP, MCF-7, T-47D and EVSA-T cell lines (at protein level).
Function	function:Protein transport. Regulator of membrane traffic from the Golgi apparatus towards the endoplasmic reticulum (ER). Has a low GTPase activity., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the small GTPase superfamily. Rab family., subunit:Interacts with SCYL1BP1. Interacts with VSP52 and RABGAP1. Isoform 1 interacts with RAB6KIFL but not isoform 2. Isoforms 1 and 2 interact with BICD1 and BICD2. Interacts with TMF1.,tissue specificity:Ubiquitous.,



Background

function:Protein transport. Regulator of membrane traffic from the Golgi apparatus towards the endoplasmic reticulum (ER). Has a low GTPase activity., PTM:Phosphorylated upon DNA damage, probably by ATM or ATR., similarity:Belongs to the small GTPase superfamily. Rab family., subunit:Interacts with SCYL1BP1. Interacts with VSP52 and RABGAP1. Isoform 1 interacts with RAB6KIFL but not isoform 2. Isoforms 1 and 2 interact with BICD1 and BICD2. Interacts with TMF1.,tissue specificity:Ubiquitous.,


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.


Products Images



Western Blot analysis of various cells using Rab 6C Polyclonal Antibody



Western blot analysis of lysates from COLO cells, using RAB6C Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).