




# RASSF7 Polyclonal Antibody

|                                    |                                                                                                                                   |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                  | YP-Ab-16232                                                                                                                       |
| <b>Isotype</b>                     | IgG                                                                                                                               |
| <b>Reactivity</b>                  | Human;Mouse                                                                                                                       |
| <b>Applications</b>                | IHC;IF;ELISA                                                                                                                      |
| <b>Gene Name</b>                   | RASSF7                                                                                                                            |
| <b>Immunogen</b>                   | The antiserum was produced against synthesized peptide derived from human RASSF7. AA range:159-208                                |
| <b>Specificity</b>                 | RASSF7 Polyclonal Antibody detects endogenous levels of RASSF7 protein.                                                           |
| <b>Formulation</b>                 | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                           |
| <b>Source</b>                      | Polyclonal, Rabbit,IgG                                                                                                            |
| <b>Purification</b>                | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.             |
| <b>Dilution</b>                    | IHC: 1/100 - 1/300. ELISA: 1/5000.. IF 1:50-200                                                                                   |
| <b>Concentration</b>               | 1 mg/ml                                                                                                                           |
| <b>Purity</b>                      | ≥90%                                                                                                                              |
| <b>Storage Stability</b>           | -20°C/1 year                                                                                                                      |
| <b>Synonyms</b>                    | C11orf13;HRAS1;HRC1;RASSF7;Ras association domain-containing protein 7; HRAS1-related cluster protein 1                           |
| <b>Observed Band</b>               |                                                                                                                                   |
| <b>Calculated Molecular Weight</b> |                                                                                                                                   |
| <b>Cell Pathway</b>                | Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Colocalizes with gamma-tubulin.                              |
| <b>Tissue Specificity</b>          | PCR rescued clones,Placenta,                                                                                                      |
| <b>Function</b>                    | similarity:Contains 1 Ras-associating domain.,                                                                                    |
| <b>Background</b>                  | similarity:Contains 1 Ras-associating domain.,                                                                                    |
| <b>matters needing attention</b>   | Avoid repeated freezing and thawing!                                                                                              |
| <b>Usage suggestions</b>           | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



## Products Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using RASSF7 Antibody. The picture on the right is blocked with the synthesized peptide.