



SERGEF Polyclonal Antibody

Catalog No	YP-Ab-16249
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA;IHC
Gene Name	SERGEF
Immunogen	The antiserum was produced against synthesized peptide derived from human SERGEF. AA range:179-228
Specificity	SERGEF Polyclonal Antibody detects endogenous levels of SERGEF protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Deafness locus associated putative guanine nucleotide exchange factor; DELGEF;Gef;Gnef;GNEFR;Guanine nucleotide exchange factor related protein; MGC128942;MGC14128;MGC14129;Secretion regulating guanine nucleotide exchange factor;TPH1;SERGEF;Secretion-regulating guanine nucleotide exchange factor;Deafness locus-associated putative guanine nucleotide exchange factor;Guanine nucleotide exchange factor-related protein
Observed Band	50kD
Calculated Molecular Weight	50kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Fetal brain,Ovary,
Function	function:Probable guanine nucleotide exchange factor (GEF), which may be involved in the secretion process.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 7 RCC1 repeats.,subunit:Interacts with SEC5. The interaction occurs only in the presence of magnesium or manganese and is stimulated by dCTP or GTP.,



Background

function: Probable guanine nucleotide exchange factor (GEF), which may be involved in the secretion process., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Contains 7 RCC1 repeats., subunit: Interacts with SEC5. The interaction occurs only in the presence of magnesium or manganese and is stimulated by dCTP or GTP.,


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.


Products Images



Western blot analysis of SERGEF Antibody. The lane on the right is blocked with the SERGEF peptide.



Western blot analysis of the lysates from HT-29 cells using SERGEF antibody.



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH 9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).