



TAGAP Polyclonal Antibody

Catalog No	YP-Ab-16253
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF;ELISA
Gene Name	TAGAP
Immunogen	The antiserum was produced against synthesized peptide derived from human TAGAP. AA range:191-240
Specificity	TAGAP Polyclonal Antibody detects endogenous levels of TAGAP protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC: 1/100 - 1/300. ELISA: 1/5000.. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TAGAP;T-cell activation Rho GTPase-activating protein;TAGAP1;FKSG15;T-cell activation GTPase-activating protein
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytosol,
Tissue Specificity	Lung,Spleen,Synovial membrane tissue,
Function	function:May function as a GTPase-activating protein and may play important roles during T-cell activation.,similarity:Contains 1 Rho-GAP domain.,
Background	This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis, celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2013],


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using TAGAP Antibody. The picture on the right is blocked with the synthesized peptide.